

























ZooMe Auto-Flex Folding Travel Scooter

Item #	Seat	Color	UOM
FLEX-AUTO	16" Folding Seat	Electric Blue	1/ea

- Fully electronic folding and unfolding in 15 seconds or less with a wireless key fob
- Aluminum alloy frame is lightweight for superior portability
- Airline safe lithium battery
- High-gloss shrouds for a protective finish
- Easy throttle control
- Delta tiller for users with limited dexterity
- Auto-folding seat with padded, adjustable armrests

- Flat-free, non-marking tires for worry-free travel
- Easy-to-adjust tiller places the scooter's controls in the ideal position
- Secondary manual-release levers and toggle switch for folding and unfolding
- Front and rear reflectors
- 2 ignition keys and 1 key fob included
- · Anti-tip wheels for additional safety

specifications			
Weight Capacity	300 lbs.	Battery Charger	2A Offboard
Top Speed	4 mph	Brakes	Electromagnetic
Estimated Range	13 miles	Freewheel Mode	Yes
Climbing Angle	6°	Front Wheels	6" x 1.5"
Turning Radius	47.2"	Rear Wheels	7" x 2.35"
Ground Clearance	1.6"	Anti-Tip Wheels	2" x .5"
Overall Dimensions	37" (L) x 17" (W) x 32.5" (H)	Base Weight	55 lbs.
Overall Dimensions Folded	17" (L) × 17" (W) × 23.5" (H)	Battery Weight (Pair)	5 lbs.
Seat Dimensions	16" (W) × 14" (D)	Heaviest Piece Weight	55 lbs.
Motor	24V x 120W	Total Weight	60 lbs.
Controller	PG S Drive 45A	Warranty on Frame	Limited Lifetime
Batteries	25.9V x 10AH Lithium	Warranty Motor Electronics Battery	24 24 12 Months