

Titan LTE Portable Powerchair

HCPCS: K0821



INSET A



INSET B




INSET C




INSET D



INSET E



INSET F



INSET G

- Disassembles for transport & storage
- Flip-back armrests with height, width and angle adjustments
- Interchangeable color panels in red & blue
- 300 lbs. weight capacity
- Flat-free, non-marking tires
- Bluetooth electronics for easy troubleshooting and diagnosing
- K0821 HCPCS (Pending PDAC Approval)

PRODUCTS	Description
TITANLTE-18FS	Titan LTE Portable Powerchair, 1/ea
SPECIFICATIONS	Description
Product Code/SKU	TITANLTE-18FS
Product Type	Power Wheelchair
Seat Width	18"
Seat Depth	18"
Seat to Floor Height	20" - 22"
Turning Radius	27"
Steering Type	Directional Joystick
Max Speed	4.2 mph
Anti Tip Wheels	Standard
Batteries	(2) 12V x 20AH
Freewheel Mode	Yes
Brakes	Electromagnetic
Casters	8" x 2"
Charger	3A Offboard
Climbing Angle	6°
Controller	Dynamic Linx 40A
Maximum Range	8 miles
Motor	200W x 3900 rpm
Drive Wheels	9" x 3"
Ground Clearance	2.5"
Heaviest Piece	39 lbs.
Base Weight (w/out battery)	70 lbs.
Seat Weight	30 lbs.
Battery Weight	22 lbs.
Base Weight	92 lbs.

SPECIFICATIONS

	Description
Warranty	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries

ACCESSORIES

	Description
SF8010	Oxygen Tank Holder, 1/ea
SF8020	Rear Basket, 1/ea
SF8030	Cane/Crutch Holder, 1/ea