



Item #	Seat	UOM
2850HD-22	22" Captain Seat	1/ea
2850HD-24	24" Captain Seat	1/ea

accessories		
TRID-2A37	Swing-away Controller Arm - Trident	1/ea
AH1000	Cane Holder Set	1/ea
AH2000	Oxygen Cylinder Holder	1/ea
AB1120	Backpack	1/ea
AZ0060	Cup Holder	1/ea
AB1010	Armrest Bag	1/ea



HCPCS
K0825



Trident HD Heavy-Duty Powerchair

- Shroudless design offers easy access to internal components and batteries.
- Semi-Reclining captain seat with adjustable headrest is height and depth adjustable.
- Armrests are padded, adjustable and removable.
- Available with 22" or 24" width Captain's seat.
- Quick-release lever allow for easy seat removal.
- Flip-up footplate includes rubber bumper to prevent damage and is height, depth and angle adjustable.
- Positioning belt included.
- Programmable Controller.
- Black, flat-free non-marking wheels.
- Large, front-anti-tip wheels feature suspension for safety, smoother ride and enhanced obstacle handling.
- Joystick is protected by metal casing.
- Easy free-wheel operation.
- Slide-out battery tray for easy access.
- Optional swing-away controller arm sold separately.
- Optional elevating legrests sold separately.

specifications			
HCPCS	K0825	Charger	5A
Weight Capacity	450 lbs.	Brakes	Electromagnetic
Top Speed	4 mph	Caster Wheels	10" x 3"
Maximum Range	15 miles	Drive Wheels	14" x 3"
Turning Radius	29"	Base Weight	98 lbs.
Climbing Angle	6°	Battery Weight (pair)	66 lbs.
Ground Clearance	2.5"	Seat Weight	48 lbs.
Dimensions	48" x 29.5"	Total Weight	212 lbs.
Seat Dimensions	22" x 20"; 24" x 22"	Warranty on Frame	Limited Lifetime
Controller	PG VR2 70A	Warranty on Electronics, Controller, Drive-Train Components	14 months
Motor	24V x 450W x 2700RPM	Warranty on Batteries	6 months
Batteries (pair)	12V x 50AH		