



Item #	Seat	UOM
VENTURA418FS	18" Folding Seat	1/cs



HCPCS
K0807



A



B



C



D

Ventura 4 Midsize Scooter

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Seat position can be adjusted forward or backward
- Front and rear suspension
- Programmable controller
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Delta tiller
- Easy free wheel operation

specifications

HCPCS	K0807	Charger	4A Offboard
Weight Capacity	400 lbs.	Brakes	Electromagnetic
Top Speed	6 mph	Rear Wheels	10" x 3" Flat-free
Maximum Range	16 miles	Base Weight	99 lbs.
Turning Radius	54"	Battery Weight (pair)	28 lbs.
Climbing Angle	7.5°	Heaviest Piece Weight	70 lbs.
Ground Clearance	2.5"	Total Weight	153 lbs.
Dimensions	48" (L) x 25" (W)	Seat Weight	26 lbs.
Floor-to-Seat Height	22" - 24"	Front Wheels	10" x 3" Flat-free
Seat Type	Folding	Warranty on Frame	Lifetime
Seat Dimensions	18" (W) x 17.5" (D)	Warranty on Electronic Controller/Drive Train Components	24 Months
Motor	24V x 350W x 5100 rpm	Warranty on Batteries	12 Months
Batteries (pair)	12V x 33AH		