

Item #	Battery	Color	UOM
ZOOME3	10AH	Red	1/ea



## ZooMe 3-Wheel Recreational Scooter

- Powerful hub driven motor providing speeds up to 15 mph.
- Folding tiller for easy transportation in a vehicle, or for storage.
- 2 keys included.
- Hand brake with locking mechanism for parking on inclines.
- Forward and reverse controls are easily changed with the switch of a button.
- Fast and slow speed settings easily changed with the switch of a button.
- Front basket for storage.
- Pneumatic tires for the smoothest ride possible.
- Extra bright LED headlight.
- Anti-tip wheel.

specifications			
Weight Capacity	300 lbs.	Battery Charger	Offboard Smart Charger
Top Speed	15 mph	Brakes	Hand Brake
Estimated Range	20 miles	Freewheel Mode	No
Maximum Climbing Angle	10°	Front Wheels	16" x 2"
Turning Radius	22"	Rear Wheels	4" x 2"
Ground Clearance	4.25"	Anti-Tip Wheels	Standard
Overall Dimensions	41" (L) x 27" (W) x 28.5" (H)	Base Weight	78 lbs.
Seat Dimensions	14.5" (W) x 10" (D)	Warranty	12 Months on Main Frame, 12 Month on Electrical Harness, Controller, Motor/Gearbox, 6 Months on Batteries
Motor	350W Brushless Hub		
Controller	350W		