

Date _____
 Sales Representative _____
 Dealer _____
 Address _____
 City _____ State/Zip _____
 Phone _____ Fax _____
 Email _____
 Reference Information _____

Quote Order
 Account Number _____
 Purchase Order Number _____
 Ship To Address _____

 Contact _____
 Phone _____

Site Details

Select one | Indoor (Standard)
 Outdoor

Power Supply (Socket) Location:
 Top Bottom Mid Landing

Floor Construction:
 Wood Concrete Tile
 Marble Granite Steel
 Other: _____

Hand Rail Removed:


Yes No

Under Floor Heat System:

Yes No

Distance to Nearest Obstruction at Top and Bottom Step:

A. _____ inches
 B. _____ inches



Number of Floors: _____

Rail Length: _____ inches

Total Number of Risers: _____

Stair Width: _____ inches

Rise _____ inches Tread _____ inches

Are Modifications planned in staircase area (door, handrail, etc)? _____



Options

- Key Switch on Lift Wide Seat
- Additional Call/Send Additional Charging Station


Track Movement Options / Configuration

- 90° Turn
- Custom Angle Turn (0-135°)
- 180° Turn
- Custom Angle Turn (135-180°)
- Double 90° Turns
- Double Custom Radius and Angle (0-135°)
- Spiral

Lift Hand: As viewed looking up



Joystick




- Left hand Right hand

Colors

- Rail Color: Beige Black
- Seat Color: Beige Black
- Seat Cushion Color: Beige Black

Seat Configuration



Standard seat height 25"

- A _____ inches
- B _____ inches
- C _____ inches

Height


- < 5'
- 5'1" to 5'6"
- 5'7" to 6'
- > 6'1"

Weight


- 150 lb
- 151 - 200 lb
- 201 - 250 lb
- 250 - 300 lb
- 300 - 350 lb

Bulkhead height from step at lowest point? _____


Top of Stairs




1. Standard Incline Finish



2. Top Overrun
Specify x _____ inches




3. 90° Top Park
Specify x _____ inches




4. 180° Top Park
Specify x _____ inches


Bottom of Stairs




1. Standard




2. Drop Nose



3. Bottom Overrun
Specify x _____ inches



4. 90° Bottom Park
Specify x _____ inches



5. 180° Bottom Park
Specify x _____ inches

Notes/Sketch

Because all of our products are custom built and painted to fit your specific needs and application, the final product is non-returnable. However, we will work with you to resolve any problems. Cancellations of orders must be done within 6 days of sign off in order to receive a full refund of the deposit and freight balances. If a cancellation is not received within 6 days then no refund on the deposit or freight balances can be given as the order will be in production.

I have read and understand the above:

Date: _____

Print Name: _____

Sign Name: _____