

Pilot

KEEP MOVING

Navigator - E604

1-800-963-7487

www.meritsusa.com



Let the Pilot by Merits 'Navigate' you to new heights with our custom built curve stairlift

Our precision camera survey kit ensures that each lift is manufactured to the exact specifications of the staircase. Each track is 100% built to order utilizing a combination of advanced 3D tube bending and skillful craftsmanship to provide a smooth, comfortable ride.

Features include:

- ◆ Synchronous folding seat and footplate
- ◆ Smooth start and stop technology
- ◆ LED on carriage to light up your staircase
- ◆ Optional attendant control

Technical Specifications

Weight Capacity	265lbs / 120Kg
Motor	950W, 24VDC
Speed	20fpm ~ 6M/min
Drive method	Rack & pinion
Incline limits	0~52 degrees
Charger	115~230V AC, 50/60Hz
Battery	24VDC re-chargable batteries
Internal arm rest width	16" - 20"
Seat to footplate height	18.8" - 20"(not adjustable)
Remote controls	2 Reprogrammable remotes included
Compliance Standards	ASME A18.1/A17.5 & EN81-40
Safety Sensors	5 Built in safety sensors as standard
Safety devices	Overtravel prevention, motor brake, key switch, retractable seat belt, armrest and swivel seat safety switches
Warranty	5 years on the motor gear box 2 years on components

