



## AEROLUX 2C

**100%**  
Carbon fiber  
frame



**Built-in dual  
battery bays**

**Caregiver  
remote  
pairing**



## KEY FEATURES

### ENHANCED RESPONSIVE CONTROL

Precise, responsive control with smooth acceleration and braking for an effortless, stable ride.



### MEMORY FOAM SEAT CUSHION



### WIDE, COMFORTABLE SEATING

Enjoy up to 19.5" of seat width with the included armrest spacers installed.



### ULTRA SLIM CROSS-FOLD

Folds to a slim 12" width for easy storage in vehicles or at home.



### SMOOTH, STABLE RIDE

Larger, reinforced wheels provide added stability and control.



### HEIGHT ADJUSTABLE FOOTREST



## SPECIFICATIONS

|                             |                                                                                               |
|-----------------------------|-----------------------------------------------------------------------------------------------|
| SKU                         | 25MNA0050109                                                                                  |
| FRAME COLOUR                | Black                                                                                         |
| NET WEIGHT                  | 33 lbs                                                                                        |
| GROSS WEIGHT                | 36 lbs                                                                                        |
| WEIGHT CAPACITY             | 300 lbs*<br>*Current FDA rating. Independently tested for 350 lbs. FDA amendment in progress. |
| DRIVING RANGE               | 20 km / 12.5 mi                                                                               |
| INCLINE CAPABILITY          | 12°                                                                                           |
| TURNING RADIUS              | 32.5"                                                                                         |
| MAXIMUM SPEED               | 6 kph / 3.7 mph                                                                               |
| MOTOR POWER                 | 2 x 200 W Brushless                                                                           |
| UNFOLDED DIMENSIONS         | 39"L x 25"W x 35"H                                                                            |
| FOLDED DIMENSIONS           | 31.5"L x 12"W x 29.5"H                                                                        |
| FRONT / REAR WHEEL DIAMETER | 8" / 12"                                                                                      |

## SEATING DIMENSIONS

|                      |                     |
|----------------------|---------------------|
| SEAT SIZE            | 18"W x 17"D x 2.5"H |
| BACKREST HEIGHT      | 16"                 |
| BETWEEN ARMRESTS     | 18.5" – 19.5"       |
| SEAT TO FLOOR HEIGHT | 20"                 |

## BATTERY

|                  |                     |
|------------------|---------------------|
| BATTERY TYPE     | Li-Ion 24 V / 10 Ah |
| BATTERY INCLUDED | One (1)             |

## POPULAR ADD-ONS



**AEROLUX 2C - FLIGHT KIT**  
SKU: 25ASG1850200



**AEROLUX 2C - CAREGIVER REMOTE**  
SKU: 25PET0150100



**AEROLUX 2C - BATTERY**  
SKU: 25PBA0150100



**AEROLUX 2C - CUSHION SETS**  
PURPLE SKU: 25ASC0150305  
RED SKU: 25ASC0150306  
BLUE SKU: 25ASC0150302



**MULTI POCKET ORGANIZER**  
SKU: 99ASG0652100



**SIDEKICK POUCH**  
SKU: 99ASG2060100