



FLIGHTWEIGHT
BY TRAVEL BUGGY



Remarkably
light at just

24 lbs



Premium
**carbon
fiber**



Caregiver
Remote Pairing
Optional Add-On



WORLD'S **LIGHTEST** POWER CHAIR

Meet the



At just 24 lbs, the F24 is the lightest power chair on the market. Featuring a 100% carbon fiber frame, ultra-slim 12" cross-fold, and travel-ready design built for comfort and stability.

KEY FEATURES

ENHANCED RIDE STABILITY

Built for confidence on every surface, larger wheels improve comfort, while a wider stance and anti-tip wheels enhance stability and control.



WIDE, SUPPORTIVE ARMRESTS



SLIM CROSS-FOLD DESIGN



CUP HOLDER



PADDED COMFORT

A cushioned seat and backrest provide enhanced comfort, while breathable mesh promotes airflow for a more enjoyable ride.



SPECIFICATIONS

MODEL SKU	26MNA0060109
FRAME COLOUR	BLACK
NET WEIGHT	24 lbs
GROSS WEIGHT	27 lbs
WEIGHT CAPACITY	265 lbs
DRIVING RANGE	12.8 km (8 mi)
INCLINE CAPABILITY	6°
TURNING RADIUS	35.4"
MAXIMUM SPEED	6 km/h (3.7 mph)
MOTOR POWER	2 x 150 W Brushless
UNFOLDED DIMENSIONS	36.2"L x 24.4"W x 33.8"H
FOLDED DIMENSIONS	29"L x 11"W x 28"H
FRONT / REAR WHEEL DIAMETER	7" / 12"

SEATING DIMENSIONS

SEAT SIZE	15.3"W x 17"D
BACKREST HEIGHT	15.5"
BETWEEN ARMRESTS	17"
SEAT TO FLOOR HEIGHT	18.5"

BATTERY

BATTERY TYPE	Li-Ion 24 V / 10 Ah
BATTERY INCLUDED	One (1)

POPULAR ADD-ONS



CAREGIVER REMOTE

SKU: 25PET0150100



EXPLORER BACKPACK

SKU: 99ASG0900100



MULTI-POCKET ORGANIZER

SKU: 99ASG0652100



SIDEKICK POUCH

SKU: 99ASG2060100



360° SWIVEL SAFETY LIGHT

SKU: 18AET04502F0



BUMPER BLANKET

SKU: 99ASG1500100