

Navigator



FORCEMECH
INNOVATE MOBILITY FREEDOM

Forbes HEALTH

5.0 ★★★★★

- Best Lightweight Electric Wheelchairs of 2024
- BEST FOR ROUGH TERRAIN

The Navigator Series - ForceMech's Flagship Foldable Powerchairs, First Introduced in 2017, Have Set the Gold Standard for Quality and Reliability while Exemplifying Mechanical Aesthetic Beauty. They Continue to be Among The Finest Portable Electric Power Wheelchairs Available on the Market Today.



- ForceMech Yellow Aesthetic Package for High Safety and Visibility
- Dual 250watt Brushless Electric Motors
- Dual Side Loading Quicksnap Forcemax¹ Lithium-Ion Batteries
- Omni Directional Joystick Control
- Electro Magnetic Breaking System
- 2-Seconds Folding and Unfolding Process
- 400lbs Weight Capacity
- 18mi² Range in a single charge
- FAA Airtravel Approved and Cruise Ready!

1: Forcemax Lithium Battery has an 20% increase in power capacity compared to other lithium battery technologies. 2: Max Range Tested with 150lbs test subject, actual distance may vary depending on user's weight.



FORCEMECH.COM 877-90-FORCE

FORCEMECH

Forcemech Navigator Spec Sheet

Unfolding Dimensions(LxWxH)	38" x 25" x 37"
Folding Dimensions(LxWxH)	25" x 14" x 31"
Ground Clearance	3.5"
Seat Cushion Size	18" x 18" x 2.75"
Turning Radius	33"
Seat Height	18"
Armrest Height	28"
Wheel Type	Polyurethane Non Flat
Front Wheel Diameter	7 x 1.8"
Rear Wheel Diameter	12.5 x 2.3"
Chair Weight	55lbs
Maximum Loading Capacity	397lbs
Battery Specification	Lithium Ion 24V 6.6Ah x2
Motor Specification	Brushless 250Watt
Controller	All Directional Joystick
Break System	Automatic Electromagnetic
Frame Material	Aircraft Aluminum Alloy
Maximum Speed	4.5MPH
Maximum Travel Range	18 Miles

The Great
Outdoors is
Calling.



Scan QR Code
for Complete
Product
Information