



# HELIX

**FREEDOM IN EVERY TURN**



**4-WHEEL PIVOT  
STEERING SYSTEM**



**FOLDS FULLY FLAT  
FOR STORAGE**



**4-WAY  
ADJUSTABLE TILLER**

## **UNSHAKABLE STABILITY. EFFORTLESS MANEUVERABILITY.**

The HELIX is the scooter that moves like a power chair. Featuring an industry-first 4-wheel pivot steering system and our unique Wide Wheel Track™ design, experience the best-in-class stability and safety where every maneuver is tight, smooth, and effortless.

# THE SCOOTER THAT MOVES LIKE A POWER CHAIR

- Featuring a 37" tight turning radius, maneuver through tight spaces and turns with the ease and agility of a power chair.
- Best-in-class stability & safety thanks to our unique Wide Wheel Track™ design.
- 4-way adjustable tiller allows you to change the height and angle to suit your needs.
- 3 backrest height levels to choose from to support good posture.
- Enhanced user comfort thanks to our fully padded, contoured seating.
- Built-in front and rear lights for added visibility at night and when reversing, so that you can ride confidently.
- Easily folds to a fully flat position for compact storage, loading, or travel.
- Built-in tow handle on the footrest makes rolling your folded scooter effortless.
- (Add-on) 4-inch thick deluxe padded seat cushion for enhanced custom comfort.
- (Add-on) Armrests with an adjustable width of 15.75" - 21.6" offer added support.



## KEY SPECIFICATIONS

<b>SCOOTER NET WEIGHT</b>	55 LBS	<b>INCLINE CAPABILITY</b>	10°	<b>UNFOLDED DIMENSIONS</b>	40"L X 24"W X 30"H
<b>SCOOTER GROSS WEIGHT</b>	65 LBS	<b>TURNING RADIUS</b>	37"	<b>SEATING SIZE</b>	15.75" W X 14.75" D X 2.5"H
<b>MAX WEIGHT CAPACITY</b>	265 LBS	<b>MOTOR POWER</b>	250 W BRUSHLESS	<b>SEAT TO DECK HEIGHT</b>	14"
<b>MAX DRIVING RANGE</b>	12.5 MI / 20 KM	<b>BATTERY TYPE</b>	LI-ION 24 V / 5.8 AH (2PCS)	<b>FRONT WHEEL SIZE</b>	8" X 2"
<b>MAX SPEED</b>	3.7 MPH / 6 KPH	<b>FOLDED DIMENSIONS</b>	35"L X 24"W X 13.7"H	<b>REAR WHEEL SIZE</b>	8" X 2.5"

## ACCESSORIZE YOUR HELIX



MULTIFUNCTIONAL BAG



SEAT BACK BAG



DELUXE PADDED SEAT CUSHION



360° CLIP-ON PHONE STAND



ADJUSTABLE ARMRESTS



CUP & PHONE HOLDER  
NOTE: ARMRESTS REQUIRED FOR ATTACHMENT