

Jazzy®

ELITE 14

The Jazzy® Elite 14 is designed for excellent indoor maneuverability and outstanding outdoor capability. The Jazzy Elite 14 features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns around corners. Large front wheels offer extra absorption over obstacles for superb climbing capabilities along with front anti-tip wheels for transitioning over grade changes and handling various terrains with ease.

*Jazzy® Elite 14
shown in Jazzy Red
with high-back seat*



Pride®
Mobility Products Corp.

 /PrideMobility

 /PrideMobilityProductsCorp

(US) 800-800-8586 • (Canada) 888-570-1113 • gojazzy.com

SPECIFICATIONS

| | |
|-----------------------------|--|
| Weight capacity | 300 lbs. |
| Maximum speed | Up to 4 mph |
| Ground clearance | 3" |
| Turning radius | 24" |
| Overall length | 35.5" |
| Overall width | 24.25" |
| Drive wheels | 14" |
| Front wheels | 3" solid anti-tips |
| Rear wheels | 6" solid casters |
| Suspension | Limited |
| Drivetrain | Two-motor, front-wheel drive |
| Braking system | Regenerative and electro-mechanical |
| Standard electronics | 50 amp, Dynamic Shark controller |
| Battery charger | Off-board, 5 amp |
| Per-charge range | Up to 9.8 miles |
| Battery requirements | (2) 12 volt, deep cycle |
| Battery size | NF-22 |
| Battery weight | 37 lbs. each |
| Base weight | 98.5 lbs. |
| Standard seat weight | 51.5 lbs. (high-back) |
| Standard seat | 18" x 18-20" and 20" x 18-20" |
| Max. seat dimensions | 20" x 20" |
| Warranty | 5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors |

COLOR CHOICES



Jazzy
Red



Jazzy
Blue

OPTIONAL SEATING

- Synergy® Seat
- Solid Seat Pan

MODELS

| | |
|-----------------------------|-------------|
| Jazzy Elite 14-2S-SS | Code: K0822 |
| Jazzy Elite 14-2S-C | Code: K0823 |

FEATURES

- Color-through shroud
- Depth adjustable, black high-back seat
- Flip-up, removable, height and width adjustable armrest
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap Belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swingaway joystick

RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Olympian
- Commander

LIFT FINDER

Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com

Pride[®]
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • gojazzy.com



/PrideMobilityProductsCorp



/PrideMobility

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.