

The Jazzy Elite HD features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front antitip wheels for handling various terrains with ease.







/PrideMobility



 $ig(m{f}ig)$ /PrideMobilityProductsCorp

SPECIFICATIONS MODELS

SPECIFICATIONS	
Weight capacity	450 lbs.*
Maximum speed	Up to 4 mph
Ground clearance	3"
Turning radius	24"
Overall length	35.5"
Overall width	24.25"
Drive wheels	14" solid
Front wheels	3" solid anti-tips
Rear wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	60 amp, Dynamic Shark controller
Battery charger	Off-board, 5 amp
Per-charge range	Up to 9.8 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	NF-22
Battery weight Base weight Standard seat weight	37 lbs. each 98.5 lbs. 51.5 lbs. (high-back)
Standard seat	20" x 18-20"
Max. seat dimensions	24" x 24"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

- · Depth adjustable, black high-back seat
- · Flip-up, removable, height and width adjustable armrest
- · Universal, fixed, in-line joystick mount
- · Flip-up height adjustable foot platform
- · Black, flat-free, non-marking tires
- 70" lap belt

OPTIONS & ACCESSORIES

Jazzy Elite HD-2HD-SS

Jazzy Elite HD-2HD-C

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Safety flag

· Lap Belt (50", 60", 80". 90")

Code: K0824

Code: K0825

- Swing-away footrest
- · Elevating leg rests
- Residual limb support
- Heel loops
- Military patches
- Swing-away joystick

RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- · Outlander Full Platform
- Olympian
- Commander

COLOR CHOICES

OPTIONAL SEATING





Jazzy Jazzy Red Blue



 Synergy[®] Seat · Solid Seat Pan



Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com



(US) 800-800-8586 • (CAN) 888-570-1113 • gojazzy.com



/PrideMobilityProductsCorp



/PrideMobility

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.