MODEL

Jazzy Select 2S-SS	K0822
Jazzy Select 2S-C	K0823
Weight Capacity	300 lbs.
Turning Radius	22.5"
Overall Width ²	23.74"
Overall Length ²	44"
Maximum Speed ¹	4.3 mph
Ground Clearance ²	0.75 at motor bracket
Front Wheels	5″ solid
Drive Wheels	10" solid
Rear Wheels	5.75" solid
Drivetrain	Two motor, mid-wheel drive®
Braking System	"Intelligent Braking" electronic regenerative, disc park brake
Suspension Type	Active-Trac®
High-back Seat Sizes	16"x16"-18", 18"x18"-20", 20"x18"-20"
Seat-to-Floor Range	21"-22"
Per Charge Range ^{2,4}	19.4 miles
Standard Electronics	40A, PG GC3
Component Weights	Base 94 lbs. <i>(heaviest piece)</i> Seat 51.4 lbs. Batteries 25.2 lbs. each
Battery Requirements ^{3,5}	(2) 12 volt, U-1 batteries
Battery Charger	3.5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

2 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

3 AGM or gel cell type required.

4 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

5 Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or -3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



(US) 800-800-8586 • (Canada) 888-570-1113 pridemobility.com | 🖸 f 🔇 🗐

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

JAZZYSELECT_12219







AMERICA'S FAVORITE POWER CHAIR gojazzy.com | 🖸 🚯 🕥 🔞



Experience the simplicity and maneuverability of the Jazzy Select[®]. Equipped with Active-Trac[®] Suspension and dual motors, this power chair ensures a smooth, comfortable ride over uneven terrain. Mid-wheel drive offers a tight turning radius, making the Jazzy Select a perfect choice both indoors and out.

- Mid-wheel drive technology and a compact size is perfect for small spaces by offering a tight turning radius of 22.5"
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Articulating rear caster beam works in conjunction with the ATX Suspension for a smooth ride over uneven surfaces
- Easy-access free-wheel lever
- · Color-through shroud keeps the shroud looking new by reducing the visual scrapes and chips

Options

Heel loops

Accessories

Cup holder

Lap belt

Flag

• Elevating leg rests

Residual limb support

Swing away joystick

Cane/crutch holder

(50", 70", 80", 90")

Swing-away footrests

Available Color



Jazzy Red

Black vinyl, depth-adjustable and limited recline standard comfort high-back with 60" lap belt

Comfort Solid Seat Pan

- 16"x16"-18"
- 18"x18"-20"

Seat Sizes

- 16"x16"-18"
- 18"x18"-20"
- 20"x18"-20"
- Rear basket Safety flag

Military patch

Oxvaen holder

- Walker holder
- Weather cover-medium
- XLR USB

Battery

Active-Trac[®] Suspension

Includes