

# Jazzy®

## SELECT ELITE

The Jazzy Select® Elite delivers high performance and maneuverability in a compact, stylish design for your active lifestyle. In-line, front-wheel drive technology gives the Jazzy Select Elite excellent stability, provides a tight turning radius around corners for optimal indoor maneuverability, climbs obstacles and grade changes with ease and handles various terrains indoors and out for an overall solid performance. The Jazzy Select Elite comes standard with Pride's exclusive Jazzy® Armor® for today's rental environment, depth-adjustable high-back seat, flip-up, removable, height and width adjustable armrest and flip-up height adjustable foot platform.



*Jazzy Select® Elite shown in Jazzy Red with high-back seat and Jazzy® Armor® (standard)*

 /PrideMobility

 /PrideMobilityProductsCorp

**Pride®**  
Mobility Products Corp.

(US) 800-800-8586 • (Canada) 888-570-1113 • [gojazzy.com](http://gojazzy.com)

## SPECIFICATIONS

<b>Weight capacity</b>	300 lbs.
<b>Maximum speed</b>	Up to 4 mph
<b>Ground clearance</b>	1.5"
<b>Turning radius</b>	24.75"
<b>Overall length</b>	34.75"
<b>Overall width</b>	22.75"
<b>Drive wheels</b>	9" solid
<b>Front wheels</b>	3" solid anti-tips
<b>Rear wheels</b>	6" solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, in-line, front-wheel drive
<b>Braking system</b>	Regenerative and electro-mechanical
<b>Standard electronics</b>	40 amp, PG GC 3 controller
<b>Battery charger</b>	Off-board, 3 amp
<b>Per-charge range</b>	Up to 15 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	U-1
<b>Battery weight</b>	24.5 lbs. each
<b>Base weight</b>	68.5 lbs.
<b>Standard seat weight</b>	43 lbs. (comfort, high-back)
<b>Standard seat</b>	18" x 18-20"
<b>Max. seat dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## COLOR CHOICES



Jazzy  
Red



Jazzy  
Blue

## OPTIONAL SEATING

- Synergy® Seat
- Solid Seat Pan

## MODELS

<b>Jazzy Select Elite-2S-SS</b>	Code: K0822
<b>Jazzy Select Elite-2S-C</b>	Code: K0823

## FEATURES

- Color-through shroud
- Jazzy® Armor® includes shroud and controller guards to protect from daily wear and tear
- Depth adjustable, black high-back seat
- Flip-up, removable, height and width adjustable armrest
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap Belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick

## RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Olympian
- Commander

**LIFT FINDER**

Use the Pride Lift Finder to find a compatible lift for your mobility product:

[prideliftsandramps.com](http://prideliftsandramps.com)

**Pride**®  
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • [gojazzy.com](http://gojazzy.com)



/PrideMobilityProductsCorp



/PrideMobility

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.