

# Jazzy<sup>®</sup>

## SPORT PORTABLE

The Jazzy<sup>®</sup> Sport Portable provides a blend of power, performance, style and portability. In-line, front-wheel drive technology provides excellent stability and maneuverability for a solid performance indoors and out while simple disassembly makes it convenient to transport.



*Jazzy<sup>®</sup> Sport Portable  
shown disassembled*

**Pride<sup>®</sup>**  
Mobility Products Corp.

 /PrideMobility

 /PrideMobilityProductsCorp

(US) 800-800-8586 • (Canada) 888-570-1113 • gojazzy.com

## SPECIFICATIONS

<b>Weight capacity</b>	300 lbs.
<b>Maximum speed</b>	Up to 4 mph
<b>Ground clearance</b>	1.5"
<b>Turning radius</b>	25.5"
<b>Overall length</b>	34.75"
<b>Overall width</b>	22.75"
<b>Drive wheels</b>	9" solid
<b>Front wheels</b>	3" solid anti-tips
<b>Rear wheels</b>	6" solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, in-line, front-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC 3 controller
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 3 amp
<b>Per-charge range</b>	Up to 11 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	U-1
<b>Battery weight</b>	24.5 lbs. each
<b>Base weight</b>	47.6 lbs.
<b>Standard seat weight</b>	43 lbs.
<b>Standard seat</b>	18" x 18-20"
<b>Max. seat dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## COLOR CHOICES



Jazzy  
Red

Jazzy  
Blue

**Pride**<sup>®</sup>  
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • [gojazzy.com](http://gojazzy.com)



/PrideMobilityProductsCorp



/PrideMobility

## FEATURES

- Easily disassembles into 5 pieces for convenient transport
- Black vinyl, depth adjustable high-back seat
- Flip-up, removable height and width adjustable armrest
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Saddlebag
- Lap Belt (50", 70", 80", 90")
- Height-adjustable footrest bracket
- Swing-away footrest
- Elevating leg rests
- Swing-away joystick

## RECOMMENDED LIFTS

- Backpacker<sup>®</sup> Plus
- Backpacker<sup>®</sup> AVP
- Outlander Full Platform
- Olympian
- Commander

**LIFT FINDER**

Use the Pride Lift Finder to find a compatible lift for your mobility product:

[prideliftsandramps.com](http://prideliftsandramps.com)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.