




Protect your investment with the Pride® Power Lift Recliner Care Kit



LC-250 shown in Cloud 9 Walnut fabric

MADE IN USA

Classic

COLLECTION

The Classic Collection from Pride® Power Lift Recliners delivers superb value and excellent comfort with features like an innovative arm design and overstuffed fiber-filled backs. In addition, Classic Collection power lift recliners feature patented, quiet and smooth operation and stylish designs to add to any home décor.

AMERICA'S #1
POWER LIFT RECLINER

POWER
RECLINERS

1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)
pridemobility.com prideliftchairs.com

Pride
Mobility Products Corp.
— LIVE YOUR BEST —

Classic COLLECTION

AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS



LC-250 shown in
Cloud 9 Pacific fabric



SPECIFICATIONS

| | Medium |
|-----------------------|--|
| | LC-250 |
| | 3-Position Chaise Lounger 5'4" - 5'10" |
| Weight Capacity | 375 lbs. |
| Seat-to-Floor | 20" |
| Seat Depth | 20" |
| Seat Width | 20" |
| Top of Back to Seat | 26" |
| Back Style | Blown Fiber Split Back |
| Fully Padded Chaise | Yes |
| Heat & Massage Option | No |
| Total Comfort Seating | Yes |
| Distance From Wall | 18" |
| Overall Chair Width | 34" |
| Chair Weight | 118 lbs. |
| Warranty | 2-Year Limited |

CLASSIC COLLECTION FABRIC SELECTION

QUICK SHIP

CLOUD 9



BLACK CHERRY



PACIFIC



STONE



WALNUT

LEXIS STA-KLEEN®

sta•kleen.
PERFORMANCE FABRIC



CHESTNUT



MUSHROOM



Look for the Quick Ship
Indicator for the shortest
lead times available.

BENEFITS

- Standard head & arm covers
- Standard large dual pockets
- Removable back
- Transformer features self-diagnostic electronics with 9' cord
- Quiet & smooth lift system
- Wooden frame composed of engineered furniture grade laminate/hardwood for added strength and durability
- Heavy-duty lift actuator & scissor mechanism
- Fire retardant (California specifications; bulletin 117-2013)
- Integrated battery backup in case of power failure¹
- Lifetime warranty on steel lift frame
- 2-year limited warranty

OPTIONS*

- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Non-skid leg levelers

Additional lead time may
apply for some options.

*Additional charges apply for all options.

¹. Batteries not included.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.