

Luggie[®]
Classic

Smart Fold, Easy Go

| | |
|----------------------------|---|
| LENGTH | 39 in (99.06 cm) |
| WIDTH | 18 in (45.72 cm) |
| WEIGHT (SCOOTER) | 52 lbs (23.59 kg) |
| WEIGHT (SCOOTER+BATTERIES) | 71.36 lbs (32.37 kg) |
| WHEELS/TIRES | Front wheels: 2"X 6" PU solid Rear wheels: 2"X 7" PU solid |
| TURNING RADIUS | 40.85 in (103.76 cm) |
| GROUND CLEARANCE | 2.5 in (6.35 cm) |
| SEAT | Foldable standard seating |
| MAX SPEED | 4 miles/hr (6.44 kph) |
| BATTERY | 2 x 12 AH AGM non-spillable batteries |
| WEIGHT CAPACITY | 250 lbs (113.4 kg) |
| DRIVE SYSTEM | 24-volt DC motor, sealed mini-transaxle |

The information contained herein is correct at the time of publication. Factors such as battery charge, battery condition, user weight, tire pressure, and terrain type will cause variation to maximum range and speed. Due to continual product improvement and manufacturing tolerances; length, width, turning radius and ground clearance can be subject to variance of +/- 3%. We reserve the right to alter specifications without prior notice.



Affordable • Compact • Simple

Easy-folding mobility scooter with tight turning radius for effortless maneuverability. Conveniently designed for daily use.

freerider[®]USA 

Toll Free: (888) 957-1888 | Se Habla Español | 我們說中文
Email: info@freeriderusa.com | Website: www.freeriderusa.com
8696 Utica Ave, Rancho Cucamonga, CA 91730