



2 0 2 6:

WHILL introduces a new lineup of ultra-lightweight mobility products



14+
Years in
Market

Sold in 20+
Countries &
Regions

Targeted
Advertising
Drives Qualified
Leads

Dealer
Support
7 days/week

RECOGNITION & CUSTOMER VOICE

~90% of users say
WHILL's design
changes perception

96% report
increased
independence



Users describe WHILL
as innovative, sleek,
and life-changing

“The Tesla of
Wheelchairs”

— Industry Media

Awards

OCES

FAST
COMPANY



TIME



EXPANDED LIGHTWEIGHT PORTFOLIO

ORDER
NOW!



MODO BREEZE

Frame:	Carbon fiber, quickly folds for transport and storage
Total weight:	28.6 lbs. (without battery and cushions)
Battery:	10Ah / 24V / 240Wh lithium-ion, safe for airline travel
Driving range:	Up to 9.5 miles per charge
Max speed:	3.7 mph
Max weight capacity:	300 lbs.
2" contoured pincore seat cushion:	18" W x 16" D

Coming
FALL
2026




WHILL MODEL C LITE

Frame:	Folding, electric wheelchair with 100% recycled carbon fiber frame
Total weight:	39.2 lbs. (without battery and cushions)
Standard battery:	10.6Ah / 25.6V / 272Wh, lithium-ion, safe for airline travel
Long range battery:	19.2Ah / 25.8V / 496Wh
Driving range:	Standard battery: up to 14.5 miles Long range battery: up to 24 miles
Max speed:	3.7 mph
Max weight capacity:	300 lbs.
Optional loading arm for extra leverage	
Bluetooth app for remote driving	

WARRANTY: Frame: 5 years/ Seat and back cushion: 2 years/ Electrical and mechanical components, batteries, and tires: 1 year. For complete warranty terms and conditions, visit whill.inc/us/warranty-policy

MARKETING & ADVERTISING

- Exclusive partner positioning & geo-targeted campaigns
- Full-funnel marketing support & lead generation
- Sales, product & technical training with certification
- Dedicated launch program: onboarding, incentives
- Unified in-store displays
- Ongoing merchandising, remarketing & dealer support



NEXT STEPS

01

Contact Adil Dean

Learn more about WHILL
adil.dean@WHILL.inc

02

Training & Onboarding

Complete certification for sales and technical staff

03

Marketing Setup

Launch geo-targeted campaigns & store finder

04

Grow your business with





modo
B R E E Z E

Ultra-Lightweight Folding Electric Wheelchair

The NEW Modo Breeze is the perfect power chair for people who need first-time or occasional mobility support.

Only
\$1,999
MAP



Ultra-light & compact



Sleek carbon fiber frame



Color LCD screen



Reinforced footrest and suspension

LED Headlight + USB Charging Port

2" Contoured Pincore Cushion (18" W x 16" D)

Calf Strap

Airline-Friendly 10Ah Li-Ion Battery

Reinforced Footrest Support

Front Suspension Springs

Rear Storage Pocket

Adjustable Back-Support Tension Straps

Adjustable Lap Belt

Rear Wheel Fenders

Anti-Tip Safety Bars

Spare Battery Bracket



HIGHLIGHTED FEATURES

- **Unit weighs 28.6 lbs.** without battery and cushions
- **300 lbs. weight capacity**
- **Reinforced footrest support** adds strength and stability for long-term durability
- **2" contoured pincore seat cushion** enhances comfort and air circulation for all-day comfort
- **Airline-friendly 10Ah Li-ion battery**
- **Programmable quick disconnect joystick** features a color LCD screen with built-in headlight + USB charging port
- **Collapsible push handles and flip-up footrest** for easy storage
- **X-frame design** allows the unit to fold nearly flat in a vertical position
- **Flip-up armrests and footplates** for easy front and side entrance

SPECIFICATIONS

Frame:	Carbon Fiber
Battery:	10Ah / 24V / 240 Wh, lithium-ion battery
Battery weight:	3.9 lbs. (1.8 kg)
Driving range:	Up to 9.5 miles
Motor:	180 W (x2); Brushless DC
Brake:	Electromagnetic Brake
Max speed:	3.7 mph
Seat size:	18" W x 16" D
Front wheel:	6.5" solid tire
Rear wheel:	11" solid tire
Size (unfolded):	38.97" x 20.86" x 34.25" (99*53*87 cm)
Size (folded):	29.92" x 8.66" x 27.16" (76*22*69 cm)
Max incline:	12 degrees



Contact Adil Dean to learn more

P: +1 (727) 303-6430

E: adil.dean@WHILL.inc

Scan QR Code to order

WARRANTY:

Frame: 5 years/ Seat and back cushion: 2 years/ Electrical and mechanical components, batteries, and tires: 1 year. For complete warranty terms and conditions, visit whill.inc/us/warranty-policy