

Outlander DE

The Outlander DE is a dual-entry, three-rail scooter lift compatible with both 3 & 4 wheel small to mid-size scooters. The lift features an integrated automatic scooter lockdown arm designed to provide a secure self-locking mechanism for compatible mobility scooters. The simple hitched-based installation requires little to no modification to your vehicle or mobility device and allows for easy removal. The drive-on/drive-off platform enables loading from either end of the lift and takes the hassle out of loading your mobility scooter in a busy parking area.

FEATURES AND BENEFITS

- · 350 lbs. lifting capacity
- · Compatible with Class II, III, or IV hitches
- LED puddle light for ease of use at night
- Toggle switch for ease of up and down use
- Integrated license plate bracket
- · Phosphate undercoating for weather resistance
- · LED indicator for On/Off
- Standard hardwire kit
- Manual override with tool kit

BASE WEIGHT ⁴	84 lbs.
MAXIMUM SCOOTER/PC OVERALL WEIGHT	350 lbs.²
MAXIMUM PRODUCT WIDTH	27"
OVERALL PLATFORM LENGTH	51", 42" usable
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	21"
MOTOR	12 Volt DC
WARRANTY	3-Year Transferable
HITCH CLASS ³	II, III, IV
OPTIONS	• Battery pack w/ charger • Swing-away adaptor • Strap kit
VEHICLE PLATFORM	Full-size van, Minivan, Truck, Sedan, SUV and Crossover
LIFTS CURRENT PRODUCTS ¹²	Victory [®] 9 3 & 4, Victory [®] 10 3 & 4, Victory [®] 10 LX w/CTS, Victory [®] Sport & Maxima 3 & 4 And equivalent competitor products Except: Pursuit [®] , Pursuit [®] XL,
	Pursuit [®] Sport 36 Volt, Sport Rider

- Vehicle dependent.
- Lifting weight does not include accessories, upgraded 2.
- Seating and/or hardware. Hitch weight carrying capacity varies by vehicle. Call for details. 3.
- 4. Does not include any accessories or hardware

To find a compatible lift for your vehicle and mobility product, use our new



pridemobility.com/lift-finder



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



liftsandramps.com • pridemobility.com 800) 800-8586 (U.S.) • (888) 570-1113 (Canada)