

Outlander XL

The Outlander XL Lift is a two-rail scooter lift specifically designed for larger 4-wheel scooters. The extra-large rails and integrated automatic scooter lockdown arm is designed to provide a secure self-locking mechanism for compatible 4-wheel scooters. The simple hitched-based installation requires little to no modification to your vehicle or mobility device and allows for easy removal. The drive-on/drive-off platform enables loading from either end of the lift and takes the hassle out of loading your mobility scooter in a busy parking area.

FEATURES AND BENEFITS

- · 350 lbs. lifting capacity
- Compatible with Class III or IV hitches
- LED puddle light for ease of use at night
- Toggle switch for ease of up and down use .
- Integrated license plate bracket
- Phosphate undercoating for weather resistance .
- LED indicator for On/Off .
- Standard hardwire kit .
- Manual override with tool kit .

BASE WEIGHT⁴	94 lbs.
MAXIMUM SCOOTER/PC OVERALL WEIGHT	350 lbs.²
MAXIMUM PRODUCT WIDTH	27"
OVERALL PLATFORM LENGTH	56.75", 49" usable
OVERALL PRODUCT WIDTH	28"
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	21"
MOTOR	Fully sealed 12 Volt DC
WARRANTY	3-Year Transferable
HITCH CLASS ³	III, IV
OPTIONS	• Battery pack w/ charger • Swing-away adaptor6 • Strap kit
VEHICLE PLATFORM	Full-size van, Minivan, Truck, SUV and Crossover
LIFTS CURRENT PRODUCTS ^{1,2}	Victory® 10 LX w/CTS, Victory® Sport, Pursuit® & Maxima 4
	And equivalent competitor products
	Except: Pursuit® XL, Pursuit® Sport 36 Volt, Sport Rider

- Vehicle dependent.
- 2.
- Lifting weight does not include accessories, upgraded seating and/or hardware. Hitch weight carrying capacity varies by vehicle.
- Call for details. Does not include any accessories or hardware

To find a compatible lift for your vehicle and mobility product, use our new



pridemobility.com/lift-finder



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice



ndramps.com • pridemobility.com 800-8586 (U.S.) • (888) 570-1113 (Canada)