

Model	SC7132
Weight Capacity	400 lbs.
Maximum Speed ²	Up to 9.3 mph
Ground Clearance ⁴	3.8" at transaxle
Turning Radius ⁴	69.3"
Overall Length ⁴	57.5"
Overall Width ⁴	29.7"
Seat-to-Ground Height	21.9" - 24.9"
Seat-to-Deck Height	17.9" - 20.9"
Front Tires	13" pneumatic
Rear Tires	13" pneumatic
Suspension	Front and rear (CTS Suspension)
Range Per Charge ^{2,3,7}	Up to 40 miles with upgradable dual lithium battery configuration
Total Weight ⁶	With standard battery: 242.6 lbs. Without batteries: 220.6 lbs.
Heaviest Piece When Disassembled ⁴	N/A (Unit does not disassemble)
Standard Seating ⁴	Type: Limited reclining high back with sliders Weight: 54.6 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 18"
Drive System	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements ⁵	25 volt, 50Ah lithium ion
Battery Charger	8 amp
Warranty	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 3 years

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC





DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

1. Pride[®] FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) +10%, -5%.
4. Tested in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Includes standard seat weight.
6. Battery weight may vary based on battery manufacturer.
7. With dual battery configuration.

0005408-070623

PURSUIT²

INSERT YOURSELF:
ANYWHERE




40 MILE
RANGE



LITHIUM-ION BATTERY
LONGER LIFE



9.3 MPH
TOP SPEED



FDA Class II
Medical Device¹

LIVE YOUR BEST[®]

pridemobility.com |    

PURSUIT²

Memory foam seat with 6" slider for adjustable support and premium, all-day comfort

Built to accompany you on all your adventures, big or small

Tiller adjustment lever

Front basket to hold whatever you need for the road ahead

Durable plastic shrouds

FDA Class II Medical Device¹

Pneumatic tires mounted on split rims for off-road thrills

ENGINEERED FOR ADVENTURE

The **Pursuit 2** proves the best adventures are yet to come as its fast cruising speed, weight capacity and pneumatic tires can tackle any road (or backroad).

FEATURES

- 400 lb. weight capacity
- 9.3 mph top cruising speed
- Up to a 40-mile range with upgradable dual lithium battery configuration
- Durable plastic shrouds
- Memory foam seat with 6" slider for adjustable support and premium, all-day comfort
- Infinitely adjustable tiller with USB charger
- Enhanced LED lighting for clear visibility
- Pneumatic tires for varied terrains and surfaces
- 3-year warranty on lithium battery
- Depth adjustable seat options up to 22" x 22"

AVAILABLE COLORS



ACCESSORIES

Only one rear accessory can be used at a time.

- 3" Single/Double Cane or Crutch Holder
- Front Basket
- Wishbone Crutch Holder
- Cane Crutch Holder
- Forearm Crutch Holder
- Walker Holder
- Cell Phone Holder
- Oxygen Tank Holder
- Cup Holder
- Saddlebag
- Lap Belt
- Weather Cover
- Safety Flag
- Rear Basket

Cup holder, rear basket and saddlebag can be ordered as a bundle.

MILITARY PATCHES



Infinitely adjustable tiller with USB charger

Enhanced LED lighting

Front and rear CTS suspension