

Enhance Mobility



TRIAXE

SPORT



The peppy TRIAXE Sport Scooter can travel at speeds up to 12mph and distance to 35 miles on a single charge. This light weight travel scooter folds flat for storage and transport.

TOUR



This lightweight, folding travel scooter features rear wheel drive, paddle acceleration and electromagnetic brakes. The aluminum frame adds strength, reduces weight and eliminates corrosion to provide years of reliability.

Construction	Aluminum Frame	Aluminum Frame
Folds Flat for Storage	Yes	Yes
Weight (without Battery)	44lbs	58 lbs
Speed	Up to 12mph	Up to 6mph
Weight Capacity	350lbs	375lbs
Battery	48V 12AH (travel battery available - 48V 6AH)	24V 12amp (12AH)
Brake	Dual Rear Disc Brakes	Electric Brake System
Length	34"	38"
Width	20 ½"	19"
Turning Radius	33"	27"
Range	Up to 35 miles	Up to 18 miles
Drive System	Front Wheel Drive Brushless Motor	Rear Wheel Drive Transaxle with Free Wheel Mode
Colors Available	Metallic Red, Light Metallic Blue, Green, Black, Silver, Purple (available soon)	Fiery Red, Bright Blue
Seat	13" standard, 18" option	13" standard, 18" option
Tires	8" Solid Front 9.5 Pneumatic Rear Tires with Anti-Flat	8" Solid Front 9.5 Pneumatic Rear Tires with Anti-Flat

Enhance Mobility

sales@enhancemobility.com

800-296-4151