

AeroFly™ Ultra-Lite Transport Chair

Packed with Features, Yet Easy to Push & Lift



AeroFly™ Ultra-Lite Transport Chair is **among the lightest transport chairs on the market!**

Weighs Only* 13 lb.

*Excludes rear wheels, and footrests. Total product weight is 17 lb. with rear wheels attached.

With a sleek, modern design and a patent-pending frame that weighs just 13 pounds*, the AeroFly™ Ultra-Lite Transport Chair is **among the lightest transport chairs on the market.** It's also competitively priced and incredibly easy to ship, store, and deliver. Ultra-lightweight yet ultra-capable, this versatile transport chair will appeal to end-users who prioritize style, comfort, safety, and convenience—at home, on the go, and everywhere in between.



Benefits & Features

Easy to Take on the Go

Offers stress-free portability your customers will love thanks to lightweight yet durable magnesium alloy construction, a fold-down backrest, quick-release rear wheels, and built-in carry handles.

Fully Equipped for Comfort & Safety

Includes attendant brakes and treaded tires for improved control, traction, and stability. Prioritize customer comfort with easy to remove seat and back overlay cushions, a standard 18" seat width, and adjustable-length footrests with integrated heel loops to help keep feet secure.

Protect Your Margins

Maximize your returns with MAP pricing to ensure you won't be undercut.

Stands Out to Sell More

Differentiate your offerings and draw customers to your showroom with the sleek and stylish patented frame design in a choice of two premium finishes for added choice.

Ultra-compact folding design ▶

Product Specs

Frame Material	Magnesium Alloy
Frame Color	Gray, Black
Rear Wheel Size	13", Quick-release
Seat Dimensions	18" (W) x 16" (D)
Seat to Floor Height	19"
Product Weight	17 lb.
Frame Weight (w/o rear wheels)	13 lb.
Weight Capacity	250 lb.
Warranty	Limited Lifetime

Item #	Description	UPC	UOM
AFT18BK-SF	AEROFly™ 18" Ultra-Lite Transport Chair, Black	822383033228	1/ea
AFT18GY-SF	AEROFly™ 18" Ultra-Lite Transport Chair, Gray	822383033235	1/ea