



Ultra-lightweight yet ultra-capable, the AeroFly™ Ultra-Lite Wheelchair is engineered for users who live life on the move. With a sleek, modern design and a frame that weighs just 13.3 lb.*, it's among the lightest wheelchairs in its class. It's competitively priced and incredibly easy to ship, store, and deliver.

Item #	Description	UPC	UOM
AF18FDABK-SF	AeroFly Ultra-Lite Wheelchair, 18", Black	822383030135	1/ea
AF18FDAGY-SF	AeroFly Ultra-Lite Wheelchair, 18", Gray	822383030142	1/ea
AF20FDABK-SF	AeroFly Ultra-Lite Wheelchair, 20", Black	822383033204	1/ea
AF20FDAGY-SF	AeroFly Ultra-Lite Wheelchair, 20", Gray	822383033211	1/ea

Product Specs	18"	20"
Frame Material	Magnesium Alloy	
Frame Weight (w/o rear wheels and anti-tippers)	13.3 lb.	14 lb.
Product Weight (includes rear wheels)	19 lb.	20 lb.
Seat Dimensions	18" (W) x 16" (D)	20" (W) x 16" (D)
Seat-to-Floor Height	19"	
Weight Capacity	250 lb.	300 lb.
Warranty	Limited Lifetime	

Accessories




Carrying Bag
Item #AFCARRYBAG



Cup Holder
Item #STDSI040S



Gray Cushions, 18"
Item #AEROFly-06-18 (Back)
Item #AEROFly-28-18 (Seat)



Orange Cushions, 18"
Item #AEROFly-03-18 (Back)
Item #AEROFly-25-18 (Seat)

Effortless Patient Adoption

Offers stress-free portability your customers will love thanks to lightweight yet durable magnesium alloy construction, fold-down backrest, quick-release rear wheels, and built-in carry handles.

Fully Equipped to Keep Customers Safer

Safety comes standard with included wheeled anti-tippers, plus a dual braking system, treaded tires for improved traction and stability, and European-style footrests with integrated heel loops to keep feet secure.

Protects Your Margins

Competitively priced against similarly equipped ultra-lightweight wheelchairs to maximize your returns with MAP pricing to ensure you won't be undercut.

Stands Out to Sell More

Differentiate your offerings and draw customers to your showroom with the sleek and stylish patent-pending frame design in a choice of two premium finishes.

*Transport weight, 18" model – excludes rear wheels, anti-tippers, and footrests. Total product weight is 19 lbs with rear wheels attached.