

# CarbonTouch™ Power Wheelchair

Travel-Ready Carbon Fiber Mobility



## Benefits & Features

### Ultra-Portable & Travel-Friendly

At just 36.4 pounds with a one-touch folding mechanism, and FAA-compliant lithium battery, the CarbonTouch is easy to carry and perfect for trains, planes, cruise ships, and more.

### Accommodates More Users

A 330-pound weight capacity and 20" (W) x 18" (D) seat accommodate more patients than competing models. Easy transition between powered and manual modes allows flexible use across care settings.

### Safety & Convenience Included

Anti-tippers, safety belt, and removable joystick with integrated USB charging port all come standard. Dual charging options allow for on-board and off board charging.

### Durable Design

Sleek carbon fiber frame offers superior durability while being significantly lighter than traditional materials.

## Meet the growing demand for travel-friendly mobility with the CarbonTouch Power Wheelchair.

With an ultra-lightweight carbon fiber frame and high weight capacity, this airline-ready power wheelchair brings premium comfort to more customers while maintaining exceptional portability for active lifestyles.



## Product Specs

Drive Range	10.5 miles
Top Speed	3.70 mph
Turning Radius	27.6"
Ground Clearance (Middle)	2.4"
Controls	Removable Joystick
Armrest Type	Flip-Back, Depth-Adjustable
Leg rest Type	Footplate
Seat Dimensions	20" (W) x 18" (D)
Seat to Floor Height	19"
Folded Dimensions	33.3" (H) X 23" (W) X 15.5" (D)

Battery Type (included)	Lithium Ion 24V 12Ah
Battery Charger Type	Off Board & Controller Port
Wheel Size	7" Front   8" Back
Base Weight	40.8 lb.   36.4 lb. without Battery
Battery Weight	4.4 lb.
Suspension	Yes
Weight Capacity	330 lb.
Warranty	1 year motor, controller and CPU 5 years frame, 1 year battery, 6 months wearable parts

Item #	Description	UPC	UOM
CFPC17FS	CarbonTouch™ Power Wheelchair	822383030319	1/ea

877.224.0946

www.drivemedical.com

99 Seaview Blvd., Port Washington, NY 11050

©2026 Medical Depot, Inc. dba Drive Medical, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive Medical are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.

CFPC17FS\_SS\_VER.A.03.26