

**1** Padded armrests and upright handle adjust to user height preferences

**2** Durable and comfortable padded seat

**3** Simple one-handed, side-to-side folding

**4** Includes removable zipper storage bag with strap

**5** Dual front fork for optimal stability and maneuverability

**6** Loop lock brakes for added safety and control

**7** Adjustable length dual hand-brake system

**8** Large 10" multi-terrain front wheels for indoor and outdoor use

**9** Height adjustable backrest

**10** Cane holder for added convenience



**NEW**

Includes sit-to-stand handles to promote safe transitioning from seated into a standing position.

The Drive Elevate Upright Walker/Rollator brings users increased balance and stability while promoting better posture and providing a more natural line of sight than traditional mobility aids. Height-adjustable upholstered forearm pads keep arms supported and level for comfortable use, while large, ten-inch wheels make it perfect for indoor and outdoor use. The **NEW** sit-to-stand handles promote safe transitioning from seated into a standing position. With a one-handed, side-to-side fold, the Elevate Upright Walker Rollator is the stand-up mobility aid that's easy to take on the go.

ITEM #	PRODUCT DESCRIPTION	UOM	UPC
UPWKRBL10	Upright Walker, Rollator, Blue	1/ea	822383015842
UPWKRGY10	Upright Walker, Rollator, Grey	1/ea	822383015859

PRODUCT SPECS	
Warranty	Limited Lifetime*
Minimum / Maximum Height	48" / 53"
Product Dimensions	25.8" (W) x 31.5" (D)
Product Weight	24.8 lbs.
Arm Pad Height	39" - 44"
Handle Extension Range	4" - 6"
Seat Dimensions	18" (W) x 10" (D) x 21" (H)
Bag Dimensions	15.5" (L) x 5.3" (W) x 9.6" (H)
Weight Capacity	300 lbs.
Color	Blue, Grey

BOX SPECS	
Carton Box Dimensions*	31" (L) x 12" (W) x 16" (H)
Carton Box Weight	27.6 lbs.



\*Refer to instruction guide or website for warranty details.