



Standard Colors:



Velvet Red



Coastal Blue

Optional Colors:



Crisp White



Desert Tan



Stone Grey



Orchid Pink



Canary Yellow

Buzzaround

3 wheel GB121B-STD

4 wheel GB124A-STD



The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers improved standard features to make your riding experience better than ever. It features a larger wraparound Delta tiller, higher ground clearance of 3.5", and an increased operating range of up to 10 miles on the standard 12 AH batteries. Upgrade with the new XR Extended Range Battery option for even more mileage! The Buzzaround XL also offers a new stylish, contoured stadium style seat for more comfort and increased seating tolerance. This scooter disassembles in seconds with just four easy steps making transporting it from one destination to another a breeze!



Features:

- 300 lb. Weight Capacity
- Improved Drivability
- Increased Operating Range
- Higher Ground Clearance
- LED Angle-adjustable Headlight
- USB Port on the Tiller
- Phone & Cup Holders
- Quick & Easy Disassembly

Upgrade with NEW XR Extended Range Battery Option!

(Packaged separately as an accessory but ships with the unit.)



LIGHTWEIGHT PORTABLE SCOOTER

Item	GB121B-STD	GB124A-STD
Description	XL	XL
Medicare Code	K0800	K0800
Weight Capacity	300 lbs.	300 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range up to ¹	10 miles/18 miles	10 miles/18 miles
Under Deck Clearance	3.5"	3.5"
Turning Radius	37.5"	45.5"
Type Batteries	2-12AH / 2-22AH	2-12AH / 2-22AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp	50 Amp

Dimensions

Length	41"	40"
Width	20"	20"
Height (ground to top of back of seat)	33.75"-35.25"	33.75"-35.25"
Ground to Top of Seat	21.75"-23.25"	21.75"-23.25"
Ground to Top of Deck	5"	5"
Deck to Top of Seat	16.5"-18"	16.5"-18"
Number of Seat Height Positions	3	3
Size of Increments	.75"	.75"
Front Axle to Rear Axle	31.5"	30.50"
Weight of Unit: (assembled)³	93 lbs.	98 lbs.
Front Half	32 lbs.	37 lbs.
Rear Half	24 lbs.	24 lbs.
Seat w/Arms	21 lbs.	21 lbs.
Batteries	16 lbs./30 lbs.	16 lbs./30 lbs.
Tires:	Solid, low profile	Solid, low profile
Front	7.5"	7.5"
Rear	8"	8"
Rear Anti-Tip Wheels	2"	2"

Standard Seat

Seat Back Height - no headrest	14.75"	14.75"
Width x Depth	17"x 14.5"	17"x 14.5"
Color	Black Vinyl	Black Vinyl

Warranty:


Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics ²	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

¹ Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

² Electronics warranty excludes batteries.



GB124A-STD 4 Wheel



Disassembles into 5 Pieces!

4 Wheel Model Shown



A battery charging port is conveniently located on the tiller under the control panel.



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com



www.facebook.com/goldentechnologies