

Nitro® Elite CF, Carbon Fiber Rollator

HCPCS: E0143 + E0156,E0143 + E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Luxurious super lightweight carbon fiber frame (Figure A)
- Frame weighs only 15 lb for powerful performance with swift and effortless mobility
- Opulent and removable zippered storage bag Unique bag attachment keeps bag securely in place when rollator is open or folded (Figure B)
- Brake cable inside frame for added safety (Figure C)
- Handle height easily adjusts with unique push button
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort (Figure D)
- Easily fold with one hand to ultra-compact size for storage Handle is located on seat (Figure E)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable
- EZ clip cane holder included

PRODUCTS	Description
RTL10266CF	Nitro Carbon Fiber Rollator, Black
SPECIFICATIONS	Description
Product Weight	15 lb (excluding bag)
Seat Dimensions	17.75" (W) x 8.5" (D)
Handle Height	32.5" - 37.5"
Weight Capacity	300 lb
Seat To Floor Height	20.75"
5 Year Limited Warranty	Brake and Handle Housing
Carton Shipping Weight	19 lb