

Nitro® Duet Rollator and Transport Chair

HCPCS: E0143 + E0156,E0143 + E0156




FIGURE A



FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Combines the features of a rollator and transport chair into one unit. The user can ambulate independently or be pushed safely by a caregiver
- The tool-free two-position contoured backrest is reversible depending on whether you are using the Nitro Duet as a rollator or transport chair
- Back support height easily adjusts with tool-free buttons
- Caster fork design and large 10" front casters enhances turning radius and rolling comfort
- Cross-brace design allows for side-to-side folding and added stability
- Brake cables routed inside the frame add safety and convenience

PRODUCTS	Description
RTL10266DT	Nitro Duet, Red, 1 Ea.
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	17.5" (W) x 13" (D)
Handle Height	33" - 37"
Weight Capacity	300 lb
Seat To Floor Height	20.5"
Casters	Back 6" Front 10"
Weight	22 lb