

Standard Fabrics:



Anchor



Bittersweet



Oxford



Port



Sandstorm

This fabric has a modern, updated look with rich color options. It features our exclusive SofTouch™ with Stain Defense coating.

* An upcharge will apply for all crossover fabrics.

Comforter PR531





The Comforter Series from Golden is all about finding the perfect size for the perfect fit. Our long-standing, Small-to-Tall size collection offers seven sizes to comfortably fit people with weight capacities up to 500 lbs. These power lift recliners feature a new back style with a tasteful seam design, our signature seat for the utmost in comfort and support, a new arm style with additional plush padding for more comfort and value, and extended arms to provide better stability while entering or exiting the chair. Enjoy luxurious new chenille fabrics to match any décor. Find your perfect fit!

Standard Features:

- Signature Seat
- Luxurious Chenille Fabrics
- Upgraded Stain Defense
- USB Charging Port



Quick Ship?

Check real-time inventory at webportalgt.com. (Access available for Golden retailers only.)



COMFORTER Series PR531



Model Number	PR531-PSA	PR531-MED	PR531-LAR	PR531-TAL	PR531-S23	PR531-M26	PR531-T28
Size	Petite/Small	Medium	Large	Tall	Small/Wide	Medium/Wide	Tall/Wide
Weight Capacity	375 LB	500 LB	500 LB				
Positioning	3	3	3	3	3	3	3
Comfort Zone(s)	1	1	1	1	1	1	1
Twilight Technology	NO	NO,	NO	NO	NO	NO	NO
Chaise Seat	YES						
Pocketed Coil Seat	YES	YES	YES	YES	NO	NO	NO
Back Type	Seam						
Removable Back Type	Bolt						
Head Rest Cover	YES						
Arm Rest Covers	YES						
Overall Width	31.5"	32"	34"	34.5"	34"	38"	40"
Overall Height	40"	41"	44"	47"	40.5"	43.5"	47.5"
Width Between Arms	21"	21"	23"	23"	23"	26"	29"
Floor to Top of Seat	19"	20"	21"	22.5"	19.5"	20.5"	22"
Seat Depth	19"	19"	20.5"	22"	20"	20.5"	22"
Seat to Top of Back	27"	27"	28.5"	30.5"	27"	29"	32"
Seat to Top of Arm	6"	5"	6"	6"	7"	7.5"	5.5"
Distance Required from Wall Reclined	16"	17"	19"	21"	16"	18"	21"
Extension Reclined from Back of Seat to End of Foot Rest	41°	41"	42"	47"	40"	42"	46"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES						
Warranty on Electrical Parts	3 YR						
Pro Rated Years 4 Thru 7	YES						
Warranty on Mechanical Labor	3 YR						
Warranty on Electrical Labor	1 YR						
Hand Control Routing	Inside Arm						
Extra Pocket Option	YES						
Heat & Massage Option	YES						
Articulating Headrest Option	NO	NO,	NO	NO	NO	NO	NO
Adjustable Lumbar Option	NO	NO,	NO	NO	NO	NO	NO
AutoFlex Option	NO	NO,	NO	NO	NO	NO	NO
Foot Rest Extension Option	YES						
Straight Lift Option	YES						
USB Charging Port	Hand Control						
Battery Backup	YES						
Place of Origin	USA						

SmartTek Diagnostic System STANDARD on all Golden models.





Comfort Zone 1 provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.





Sizes may vary depending on fabric, filling material,

upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.




















Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374 Fax: 800-628-5165 www.goldentech.com

