

M A X I C O M F O R T <sup>TM</sup> S E R I E S



Day Dreamer  
PR630  
shown in Coffee Bean


**Pocketed Coil  
Spring Seat  
Standard!\***



Hazelnut



Calypso



Shiraz



Coffee Bean\*  
\* upcharge applies



The DayDreamer with the exclusive PowerPillow™ provides the ultimate in position control. Recline back and adjust your head and neck using the exclusive powered head pillow. MaxiComfort combined with the PowerPillow offers maximum comfort options.





Day Dreamer



Power Articulation Head Forward



Power Articulation Head



Power Articulation Head Back

|  |              |
|--|--------------|
| Model Number   | PR630        |
| Size   | Medium       |
| Weight Capacity  | 375 Lb.      |
| Recline Positions  | MaxiComfort  |
| Standard Fabrics   | 4            |
| Fabric Options   | ALL          |
| Chaise Seat  | YES          |
| Coil Springs   | YES          |
| Back Type  | Seam         |
| Head Rest Cover  | NO           |
| Arm Covers   | NO           |
| Battery Back Up  | YES          |
| Overall Width  | 38"          |
| Overall Height   | 44"          |
| Width Between Arms   | 21"          |
| Floor to Top of Seat   | 20"          |
| Seat Depth   | 19.5"        |
| Seat to Top of Back  | 28.5"        |
| Seat to Top of Arm   | 5.5"         |
| Seat Angle Raised  | 16°          |
| Floor to Top of Foot Rest Reclined                             | 23"          |
| Floor to Top of Back Reclined                                  | 23.5"        |
| Back Angle Closed  | 105°         |
| Back Angle Reclined  | 170°         |
| Foot Rest Angle Closed   | 90°          |
| Foot Rest Angle Opened   | 175°         |
| Distance Required from Wall Reclined                           | 25.5         |
| Extension reclined from Back of Seat to end of Foot Rest       | 42"          |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | YES          |
| Warranty on Electrical Parts                                   | 3 Yrs.       |
| Pro Rated Years 4 Thru 7                                       | YES          |
| Warranty on Mechanical Labor                                   | 3 Yrs.       |
| Warranty on Electrical Labor                                   | 1 Yr.        |
| Removable Back Option Available                                | NO           |
| Smart Tek System Standard                                      | YES          |
| Customizable   | NO           |
| Left Hand Control Option Available                             | NO           |
| Extra Pocket Option Available                                  | NO           |
| Heat & Massage Option Available                                | YES          |
| Foot Rest Extension Option Available                           | NO           |
| Straight Lift Option Available                                 | NO           |
| User Height (suggested)  | 5'4" - 5'10" |

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



**Battery Backup System**

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



**Smart Tek**

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68  
Listing # E112417  
FNAE 80-214A

