

M A X I C O M F O R T TM S E R I E S



Relaxer
PR756L



Relaxer
PR756MC



PR756L
FABRICS



Hazelnut



Sterling



Copper



Calypso



Shiraz

PR756MC
FABRICS



Hazelnut



Sterling



Copper



Calypso



Shiraz



Buckskin*



Coffee Bean*

* upcharge applies

Relaxer
PR756MC
PR756L

The Relaxer is the original powered lift chair recliner featuring MaxiComfort® Ultimate Recline Technology standard! The Relaxer model is a true Golden customer favorite! Featuring the most luxurious four-pillow waterfall back design, and overstuffed, super plush rounded armrests, the Relaxer is available in two sizes – Medium and Large – for the absolute perfect fit. Every MaxiComfort chair comes standard with our exclusive AutoDrive™ 3.0 hand control to make reclining the chair into your favorite position a breeze!



Relaxer



Relaxer
PR-756MC **Relaxer**
PR-756L

Model Number	PR-756MC	PR-756L
Size	Medium	Large
Weight Capacity	375 LB	375 LB
Recline Positions	MaxiComfort	MaxiComfort
Standard Fabrics	7	5
Fabric Options	ALL	ALL
Chaise Seat	YES	YES
Coil Springs	NO	NO
Back Type	4 Pillow	4 Pillow
Head Rest Cover	YES	YES
Arm Covers	NO	NO
Battery Back Up	YES	YES
Overall Width	36.5"	39"
Overall Height	44"	46"
Width Between Arms	20"	22"
Floor to Top of Seat	19"	20.5"
Seat Depth	21"	22"
Seat to Top of Back	29"	30"
Seat to Top of Arm	8"	8"
Seat Angle Raised	15°	15°
Floor to Top of Foot Rest Reclined	23"	24"
Floor to Top of Back Reclined	21"	22"
Back Angle Closed	100°	100°
Back Angle Reclined	174°	174°
Foot Rest Angle Closed	90°	90°
Foot Rest Angle Opened	180°	180°
Distance Required from Wall Reclined	25"	27"
Extension Reclined from Back of Seat to End of Foot Rest	43"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.
Removable Back Option Available	NO	NO
Smart Tek System Standard	YES	YES
Customizable	NO	NO
Left Hand Control Available	YES	YES
Extra Pocket Option Available	YES	YES
Heat and Massage Option Available	YES	YES
Foot Rest Extension Option Available	YES	YES
Straight Lift Option Available	NO	NO
User Height (suggested)	5'4"-5'10"	5'11"-6'2"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68
Listing # E112417
FNAE 80-214A

