

# FEATURES AND SPECIFICATIONS

Specifications	Details
Scooter Size Setup	39" H x 38" L x 22" W
Scooter Size Folded with Seat	30.5" H x 38" L x 22" W
Scooter Size Folded without Seat	17" H x 38" L x 22" W
Scooter Weight Complete	39.5lbs Total 27 lbs Heaviest Part
Seat Weight	8.5 lbs
Battery Weight	4 lbs
Maximum Load Weight	300 lbs loaded or 136 kg
Scooter Speeds (3)	3, 5 and 7 mph
Ground Clearance	9" front - 3" rear
Battery	36V, 8 amp hr 288 watt Li-Ion Rechargeable
Battery Certification	UN PT. III 38.3 for Airline Travel
Battery Charger	Input Voltage : AC100V-240V
Battery Charging Time	~1 hour (~5 hours for initial charge)
Frame Materials	Stainless Steel
Brakes	Front Disk Brake
Drive	Front Wheel Drive

